

EMERGENCY RESPONSE KIT CABINET

MACHLAB[®] Emergency Response Cabinet have distinctive features that are designed to provide safe storage of emergency response kits to meet the health and safety requirements of the laboratory today. All laboratories shall be equipped with sufficient safety kits that are easily accessible to response to any emergency promptly.



CONTENTS

First Aid Kit (MS-S)

- 1 x Acriflavine solution – 30ml
- 1 x Liniment methyl salicylate – 30ml
- 1 x Antiseptic cream – 15gm
- 1 x Antiseptic lotion – 60ml
- 1 x Absorbent cotton wool (w/box) – 5gm
- 1 x Gauze swabs (w/box) – 8pcs
- 1 x W.O.W bandage 5cmx5.5cm
- 1 x Medicated oil – 3ml
- 10 x Plaster strips (fabric)
- 10 x Paracetamol tablets (UPHAMOL)
- 2 x Sterile wipe
- 1 x Adhesive plaster roll
- 1 x Scissor non-foldable
- 1 x Tweezer
- 1 x First aid book
- 1 x Usage card

Chemical Spill Kit (SK571717)

- 1 x 20 litre pail
- 2 x Chemical sorbent sock (1pc able to absorb >8.25 litres)
- 2 x Chemical sorbent pillow (1pc able to absorb >2.5 litres)
- 5 x Chemical sorbent pads (1pc able to absorb >0.725 litres)
- 2 x Hazardous material disposable bags 24" x 30"
- 2 x Cable ties
- 1 x Anti-fog goggle
- 1 x Nitrile gloves

Fire Blanket (FB44CD)

- 1 x Fire blanket (4' x 4')

Emergency Response Cabinet

- Tempered safety glass window
- Prominently labelled
- Wall-mounted at eye level, away from any obstruction

FEATURES

- Durable **lab-grade construction**
- **100 microns** powder coating finish
- **Corrosion resistance** certified
- **Safety Yellow** colour for emergency access
- **Prominent safety label** and brightly coloured

SPECIFICATIONS

- **Certified quality** (ISO2178, ISO2813, ISO2409, ASTM 117, ASTM D3363)
- Manufactured of **heavy-duty steel** (24-gauge, cold-rolled)
- Strong **adjustable shelf**
- Stainless steel/Aluminium alloy fittings (**Antirust**)

MODELS

	Size (mm)			Packing Size (mm)	Packing Volume (m ³)
	W	D	H		
ERK-4575	450	400	750	470 x 420 x 770	0.152

