



## HAZARDOUS STORAGE CABINET

MACHLAB<sup>®</sup> safety cabinets have distinctive features that meet international regulatory and compliance requirements for safe chemical storage and protection.



Standard

Under-bench



### FEATURES

- **High performance design** to meet OSHA and NFPA standards for safe storage of hazardous materials.
- **Maximum protection and security** with fail-safe 3-point self-latching lock and full continuous piano hinged doors.
- Constructed of 1.0-1.2mm thick double walled, fully welded steel with 38mm of insulating air space for **fire resistance**.
- **Prominent warning label** in multiple languages.

### COMPLIANCE & REGULATORY REQUIREMENTS

- **OSHA** 29 CFR 1910.106.
- **NFPA** Code 30 Section 5.9.3.
- Federal, state, and local regulations may apply.
- Contact your local Fire Marshall for guidance.

### MODELS

	Capacity		Size (mm)			Adjustable Shelves	Packing Size (mm)	Packing Volume (m <sup>3</sup> )
	Gallons	Litres	W	D	H			
<b>SCH22</b>	22	83.0	889	559	889	1	909 x 579 x 909	0.478
<b>SCH30</b>	30	114.0	1090	450	1120	1	1110 x 470 x 1140	0.595
<b>SCH45</b>	45	170.0	1090	450	1650	2	1110 x 470 x 1670	0.871
<b>SCH60</b>	60	227.0	850	850	1650	2	870 x 870 x 1670	1.264
<b>SCH90</b>	90	341.0	1090	850	1650	2	1110 x 870 x 1670	1.613