


REAGENT STORAGE CABINET – RHC-3518

MACHLAB[®] Reagent storage cabinets have distinctive features that are designed to provide convenient storage for your chemical reagents and meet international laboratory safety standards.



FEATURES

- Durable **lab-grade construction**
- **100 microns** powder coating finish
- **Corrosion resistance** certified
- **Fail-safe heavy-duty** sliding mechanism
- **Exhaust ventilation collar** add-on available

FINISH OPTIONS

Powder coated steel components or stainless-steel finish. Custom colour matching available upon request.



SPECIFICATIONS

- **Certified quality** (ISO2178, ISO2813, ISO2409, ASTM 117, ASTM D3363)
- Manufactured of **heavy-duty steel** (22-gauge, electro-galvanised)
- Stainless steel/Aluminium alloy fittings (**Antirust**)

MODELS

	Size (mm)			Packing Size (mm)	Packing Volume (m ³)
	W	D	H		
RHC-3518	350	600	1800	370 x 620 x 1820	0.418