

WALL MOUNTED SHELVING 2-tier



Simple. Durable. Functional. Machlab's quality is reflected in the seamless yet functional design; durable yet versatile construction; that is both convenient to use and adjustable to fit various reagent bottle sizes.

FEATURES

- Durable **lab-grade construction**
- **100 microns** powder coating finish
- **Corrosion resistance** certified
- **Safety** guardrail available on request
- Best industry **3 years warranty**

SPECIFICATIONS

- **Certified quality** (ISO2178, ISO2813, ISO2409, ASTM 117, ASTM D3363)
- Manufactured of **heavy-duty steel** (16 & 22-gauge)



	750mmW	900mmW	1050mmW	1200mmW
2-tiers	RS2-750	RA2-900	RA2-1050	RA2-1200
3-tiers	RS3-750	RA3-900	RA3-1050	RA3-1200