


GLASS SWING DOOR HIGH CABINET HCG

The Machlab design team has recreated the high cabinet to meet international ISO standards of quality and durability. The cabinet features solid construction with reinforced shelves, doors and carcass. All fittings are made of stainless steel and aluminium alloy to withstand harsh working environments. Functional yet aesthetically pleasing – it exuberates a smart, clean and modern look.



FEATURES

- Durable **lab-grade construction**
- **100 microns** powder coating finish
- **Corrosion resistance** certified
- **3 point security** locking system
- Extra adjustable shelves available
- Optional **PP spill tray** available
- Best industry **warranty (3 years)**

FINISH OPTIONS

Powder coated steel components or stainless steel finish.
Custom colour matching available upon request.



SPECIFICATIONS

- **Certified quality** (ISO2178, ISO2813, ISO2409, ASTM 117, ASTM D3363)
- Manufactured of **heavy-duty steel** (22-gauge, electro-galvanised)
- U-channel **reinforced adjustable shelves**
- Stainless steel/Aluminium alloy fittings (**Antirust**)
- Framed, 5.0mm **tempered safety glass**

MODELS

	W	D	H
HCG-8045	800	450	1800
HCG-9045	900	450	1800
HCG-1045	1000	450	1800