



## POLYPROPYLENE (PP) SAFETY STORAGE CABINET

MACHLAB® safety cabinets have distinctive features that meet international regulatory and compliance requirements for safe chemical storage and protection.



### FEATURES

- **Leak-proof** adjustable shelves.
- Padlock hasp for **security protection**.
- 2" liquid tight **containment sump** to contain any accidental leaks and spills internally.
- **Prominent warning label** in multiple languages.
- Fully PP hot weld construction for **maximum corrosion resistance**.

### COMPLIANCE & REGULATORY REQUIREMENTS

- Meets **OSHA** secondary containment requirements
- Federal, state, and local regulations may apply.
- Contact your local Fire Marshall for guidance.

### MODELS

	Capacity		Size (mm)			Adjustable Shelves	Packing Size (mm)	Packing Volume (m <sup>3</sup> )
	Gallons	Litres	W	D	H			
<b>SCP30</b>	30	114.0	1090	460	1120	1	1110 x 480 x 1140	0.607
<b>SCP45</b>	45	170.0	1090	460	1650	2	1110 x 480 x 1670	0.890
<b>SCP60</b>	60	227.0	860	860	1650	2	880 x 880 x 1670	1.293
<b>SCP90</b>	90	341.0	1090	860	1650	2	1110 x 880 x 1670	1.631